


GEBRAUCHSANWEISUNG

Personal Protection



Chirurgische Gesichtsmasken mit Ohrschlaufen

So tragen Sie Ihre Crosstex-Gesichtsmaske richtig:



Step 4 for SecureFit™ Masks only

Ref.#'s GCLBL, GCLPK, GCLBL*, GCPK*, GCAPTPK*, GCIBLSF, GCILVSF, GCIPKSF, GCIBL, GCIGR, GCILV, GCIPK, GCISA, GCITE, GCITZ, GCITQ, GCIBL100, GCIBL300, GCICXB, GCICXS, GCICXZ, GCICXT, GCIPWB, GCIPWS, GCIPWZ, GCIPWT, GCICXBSF, GCIPWBSF, GPLUSBL, GPLUSKA, GPLUSPK, GPLUSWHSF, GCPBLSF, GCPLVSF, GCPPKSF, GCPBL, GCPLV, GCPPK, GCPYE, GCPBL100, GCABL*, GCAPK*, GCFCXSSF, GCFCXSFSF, GCPWSSF, GCFCXUSF, GCFCXS, GCFCXSf, GCPWS, GCFCXU, GCFCX, GCPW

Verwendungszweck:

Crosstex™ OP-Masken dienen der Infektionskontrolle, um die Kontamination durch eingeatmete und ausgeatmete Mikroorganismen zu minimieren und den Kontakt des Trägers mit Blut und Körperflüssen zu verringern. Es handelt sich um ein Einwegprodukt, das unsteril geliefert wird.

Anwendungsgebiete: Crosstex™ OP-Masken dienen der Infektionskontrolle, um die Kontamination durch eingeatmete und ausgeatmete Mikroorganismen zu minimieren und den Kontakt des Trägers mit Blut und Körperflüssen zu verringern. Es handelt sich um ein Einwegprodukt, das unsteril geliefert wird.

Kontraindikationen:

Nicht wiederverwenden, dies ist ein Einmalprodukt.

Nicht in der Nähe von entzündbaren Narkosegasen verwenden.

Einhaltung:

US-FDA 510K cleared class II medical device Surgical Mask.

US-FDA Regulation: 21 cfr 878.4040.

Meets fluid barrier protection per ASTM F2100 standards.

Meets applicable ISO 10993 Biocompatibility testing requirements.

Meets Class 1 Flammability Spread per 16 cfr 1610.

Devices are manufactured in compliance with 21 cfr Part 820 and ISO 1385:2016 quality systems.



Crosstex International, Inc.
6789 W. Henrietta Road, Rush, NY 14543 USA
585.359.0130 | 631.582.6777 | crosstex.com



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Do not reuse.



Not made with natural rubber latex.

Made in USA unless otherwise noted. *Made outside of the USA.
EARLOOPMASK-IFU / REV D - 07/2021

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
©2021 Crosstex International, Inc. All rights reserved.

CROSSTEX

Anwendbare Symbole und Definitionen Verwendung:

 crosstex.com	Consult electronic instructions for use
	Catalog number
	Single use only Do not re-use
	Non-sterile
	Manufacturer
	Unique Device Identifier
	Manufacturing lot number
	Date of manufacture
	Medical Device
	Approved for European Union, notified body
	European Union Representative